SAFECARE *One Step Rapid Test*

COVID-19 Antigen Rapid Test Kit(Saliva)

Package Insert

For Medical professional in vitro diagnostic use only.

INTENDED USE

The COVID-19 Antigen Rapid Test Kit (Saliva) is a lateral flow immunoassay intended for qualitative detection of nucleocapsid protein antigen in oropharyngeal saliva specimens from individuals suspected of COVID-19 by their healthcare provider.

The COVID-19 Antigen Rapid Test Kit (Saliva) does not differentiate between SARS-CoV and SARS-CoV-2. Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in oropharyngeal saliva during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results from patients with symptom should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

The COVID-19 Antigen Rapid Test Kit (Saliva) is intended for use by medical professionals or trained operators who are proficient in performing rapid lateral flow tests.

PRINCIPLE

The COVID-19 Antigen Rapid Test Kit (Saliva) is an immunochromatographic membrane assay that uses highly sensitive antibodies to detect SARS-CoV-2 nucleocapsid protein from oropharyngeal saliva specimens.

SARS-CoV-2 specific antibodies are immobilized onto the test region of the membrane and combined with other reagents/pads to construct a test strip.

During testing, the specimen reacts with anti-COVID-19 antibodies conjugated to colored particles and precoated onto the sample pad of the test. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the reagents in test line region. Therefore, if the specimen contains COVID-19 antigen, a colored line will appear in test line. If the specimen does not contain COVID-19 antigen, no colored line will appear in the test line regions, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

KIT COMPONENTS

Individually packed test devices	Each device contains a strip with colored conjugates and reactive reagents pre-spread at the corresponding regions
Extraction tube	For specimens preparation
Extraction buffer	For specimens extraction
Workstation	For placing extraction tube
Funnel	For specimens collection
Package insert	For operation instruction

MATERIALS REQUIRED BUT NOT PROVIDED

Timer

WARNINGS AND PRECAUTIONS


- Do not use after expiration date. Do not use if pouch is damaged or open. Do not reuse the tests.
- Do not mix components from different kit lots. Avoid cross-contamination of specimens by using a new specimens collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Humidity and temperature can adversely affect results.
- The extraction buffer contains a salt solution if the solution contacts the skin or eyes, flush with copious amounts of water.
- Discard the using testing materials in accordance with local regulations.

REAGENT PREPARATION AND STORAGE INSTRUCTIONS

Store unused test devices unopened at 4°C-30°C. If stored at 4°C-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

SPECIMEN COLLECTION AND HANDLING

- Inadequate specimen collection or improper specimen handling may yield a false result.
- Posterior oropharyngeal saliva : Perform hand hygiene with soap and water/alcohol-based hand rub. Make a "Kruuuua" noise from the throat to clear the saliva from deep throat, then spit saliva into the funnel to the fill line(300µL) on the extraction tube. Avoid any saliva contamination of the outer surface of the tube
- Squeeze the buffer vial and add all the extraction buffer into the extraction tube.
- Insert a dropper tip into the specimen extraction tube tightly. Shake and mix completely.

- Optimal timing of specimen collection: After getting up and before brushing teeth, eating or drinking.
- The test should be performed immediately after the sample is collected. Do not leave the sample at room temperature for more than 2 hours. Specimens may be stored at -20°C for up to 1 month prior to test.
- If Specimens are to be transported, they should be packed in compliance with local regulations covering the transportation of etiological agents.
- If samples are stored at -20°C they must be returned to room temperature, thawed completely and fully mixed prior to testing. The samples can be frozen and thawed once and repeated freezing and thawing should be avoided.

ASSAY PROCEDURE

Allow the test device and specimens to equilibrate to room temperature(15-30°C or 59-86°F) prior to testing. For best results, the test should be performed in one hour.

- Remove the test device from the sealed pouch.
- Reverse the specimen extraction tube, holding the specimen extraction tube upright, transfer 3 drops to the specimen well(S) of the test device, then start the timer.
- Interpret the results at 15 minutes. Do not interpret the results after 20 minutes.

INTERPRETATION OF ASSAY RESULT

POSITIVE RESULT:

A colored line appears in the control line region (C) and a colored line appears in test line region(T).

***NOTE:** The intensity of the color in the test line region will vary depending on the concentration of COVID-19 antigen in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE RESULT:

A colored line appears in the control line region (C) and no line appears in test line region(T).

INVALID RESULT:

No line appears in control line region(C). Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the procedure with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

- Internal Control:** This test contains a built-in control feature, the C band. The C line develops after adding sample solution. Otherwise, review the whole procedure and repeat test with a new device.
- External Control:** Good Laboratory Practice recommends using the external controls, positive and negative to assure the proper performing of the assay.

PERFORMANCE CHARACTERISTICS

1. Clinical study: A side-by-side comparison was conducted using the research reagent and referencing reagent. The comparative PCR sample were obtained from oropharyngeal swab, and the sample collection for the antigen test is carried out according to this instruction. Compare with RT-PCR:

		PCR result		Total
		Positive	Negative	
Safecare Test	Positive	131	1	132
	Negative	2	182	184
Total		133	183	316
Relative Sensitivity		98.50% (94.67% ~ 99.82%)		
Relative Specificity		99.45% (96.99% ~ 99.99%)		
Overall Agreement		99.05% (97.25% ~ 99.80%)		

2. Cross-reactivity: Cross-reactivity studies are performed to demonstrate that the test does not react with the following microorganisms in the table below.

Cross-Reactant	Concentration	Cross-Reactant	Concentration
Adenovirus	1.0 x 10 ⁵ TCID ₅₀ /mL	Influenza A	1.0 x 10 ⁵ TCID ₅₀ /mL
Human metapneumovirus (hMPV)	1.0 x 10 ⁵ TCID ₅₀ /mL	Influenza B	1.0 x 10 ⁵ TCID ₅₀ /mL
Rhinovirus	1.0 x 10 ⁵ TCID ₅₀ /mL	MERS	1.0 x 10 ⁵ TCID ₅₀ /mL
Enterovirus	1.0 x 10 ⁵ TCID ₅₀ /mL	Chlamydia pneumoniae	1.0 x 10 ⁵ CFU /mL
Human coronavirus OC43	1.0 x 10 ⁵ TCID ₅₀ /mL	Bordetella pertussis	1.0 x 10 ⁵ CFU /mL
Human coronavirus 229E	1.0 x 10 ⁵ TCID ₅₀ /mL	Haemophilus influenzae	1.0 x 10 ⁵ CFU /mL
Human coronavirus NL63	1.0 x 10 ⁵ TCID ₅₀ /mL	Legionella pneumophila	1.0 x 10 ⁵ CFU /mL
Human coronavirus HKU1	1.0 x 10 ⁵ TCID ₅₀ /mL	Mycoplasma pneumoniae	1.0 x 10 ⁵ CFU /mL
Human parainfluenza virus 1	1.0 x 10 ⁵ TCID ₅₀ /mL	Streptococcus pneumoniae	1.0 x 10 ⁵ CFU /mL
Human parainfluenza virus 2	1.0 x 10 ⁵ TCID ₅₀ /mL	Streptococcus pyogenes	1.0 x 10 ⁵ CFU /mL
Human parainfluenza virus 3	1.0 x 10 ⁵ TCID ₅₀ /mL	Mycobacterium tuberculosis	1.0 x 10 ⁵ CFU /mL
Human parainfluenza virus 4	1.0 x 10 ⁵ TCID ₅₀ /mL	Staphylococcus aureus	1.0 x 10 ⁵ CFU /mL
Respiratory Syncytial Virus	1.0 x 10 ⁵ TCID ₅₀ /mL	Candida albicans	1.0 x 10 ⁵ CFU /mL

3. Interference: The following endogenous interference substances were evaluated at the concentrations listed and no effect was found.

Whole blood(2%), three OTC nasal sprays(10%), three OTC nasal drop(25%), three nasal mouthwashes(25%), 4-Acetamidophenol(10mg/ml), Acetylsalicylic acid(20mg/ml), Chlorpheniramine (5 mg/ml), Dexamethorphan (10mg/ml), Diphenhydramine(5mg/ml), Ephedrine(20mg/ml), Guaiaacol glyceryl ether(20mg/ml), Oxymetazoline(10mg/ml), Phenylephrine(100mg/ml), Phenylpropanolamine(20mg/ml)

LIMITATIONS OF TEST

1.For professional in vitro diagnostic use and should be only for the detection of SARS-CoV-2 antigen, not for any other viruses or pathogens.

2.This test detects both viable (live) and non-viable SARS-CoV-2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.

3.The performance was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.

4.A negative test result may occur if the level of antigen in a sample is below the detection limit of the test. If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. As with all diagnostic tests, a confirmed diagnosis should only be made by physician after all clinical and laboratory findings have been evaluated.

5.Negative test results are not intended to rule in other non-SARS viral or bacterial infections.

6.Positive and negative predictive values are highly dependent on prevalence. False positive test results are more likely during periods of low COVID activity when prevalence is moderate to low.

7.False negative results may occur if a specimen is improperly collected, transported, or handled.





8.Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between children and adults.

9.If the differentiation of specific SARS viruses and strains is needed, additional testing, in consultation with state or local public health departments, is required.

REFERENCES

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011;81:85-164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology. 6th ed. Lippincott Williams & Wilkins, 2013:825-58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016;24:490-502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019;17:181-192.
5. Wong G, Liu W, Liu Y, Zhou B, Bi Y, Gao GF. MERS, SARS, and Ebola: the role of super-spreaders in infectious disease. Cell Host Microbe 2015;18:398-401.

INDEX OF SYMBOLS

	Do not reuse		For in vitro diagnostic use only
	Stored between 4-30°C		Consult instruction for use
	Keep away from sunlight		Lot number
	Use by date		Contains sufficient for <n> tests
	Manufacturer		Date of Manufacture
	Sterilized using ethylene oxide		Sterilized using irradiation
	Authorized Representative in the European Community		

Safecare Biotech (Hangzhou) Co., Ltd.
Building 2/203, No.18 Haishu Rd, Cangqian Sub-district
Yuhang District Hangzhou, 311121, China

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Safecare Biotech(Hangzhou) Co., Ltd.
NIC GmbH